



The **Single-Operator** Revolution

Mastering polypectomy with Holdfinger

The *Revolution* is here

Precision, control, and reliability in your hands.

Digestive endoscopy is undergoing a major evolution, driven by **breakthrough technologies** that enhance both safety and effectiveness in clinical practice. Among these innovations, Holdfinger is emerging as a game-changing device designed to transform the execution of complex procedures such as **polypectomy** and **endoscopic mucosal resection (EMR)**.

Holdfinger is an **innovative and revolutionary medical device** designed to empower endoscopists with greater autonomy and precision during procedures. This **ground-breaking tool** enables physicians to simultaneously control and maneuver endoscopic accessories while directly managing their operation, all in real time.

Traditionally, these interventions required **multiple operators** to ensure precision and stability. Holdfinger challenges this paradigm by enabling a single operator to perform advanced techniques with **greater accuracy, stability, and efficiency**.



Holdfinger changes this dynamic by placing direct control back into the hands of the endoscopist, reducing dependence on support staff and improving the overall efficiency, precision, and safety of the procedure.

01



LOCKING PIN

Holdfinger features a specific locking system, which is reversible, fast, easy and user-friendly.

02



"T" SHAPED BASEPIN

Holdfinger allows you both moving the accessory backwards and forwards and rotate the snare.

03



ADJUSTABLE BAND

The adjustable band allows to match the measure of the ring on your finger.



By using Holdfinger, the endoscopist can manage the procedure independently and use the necessary accessories during the procedure (snare, clip, etc.).

Quick first use


How to set up Holdfinger before using it




Modify the size using the adjustable band.



Pass the catheter through the slot, ensuring that the adjustable band is facing the operator.



Align the Holdfinger facing upward and position it 2 cm from the valve.



Insert your 4th finger into the Holdfinger



With Holdfinger you can rotate, move back and forth, open and close the snare inside the channel completely by yourself.



Watch the full video on YouTube



Scan the
QR Code

Or click here:

https://youtube.com/lllCq3Zg25l?si=_89__B_fUGxj73-T



Holdfinger advantages

1

Lower direct intervention costs

Reduce the financial burden of procedures, saving on both equipment and operational costs.

2

Shorter assistant time

Decrease the duration of procedures, leading to faster patient throughput and less discomfort.

3

Shorter assistant time

Decreases the duration assistants are required for each procedure, increasing the overall clinical workflow efficiency.

4

Human resources reallocation

Free up medical staff for other tasks, optimizing workforce efficiency and resource use.

4

More control for the endoscopist

Improve procedural accuracy and flexibility, enhancing outcomes through superior device maneuverability and feedback.

6

Minimized complications

Lower the risk of procedural errors and adverse events, improving patient safety and recovery time.



The Single Operator Procedure was not created to replace the role of the nurse, but to optimize the entire endoscopic workflow.

Get To Know Us Better

Do you want more info about Holdfinger and receive a sample?

Contact us at vendite@meditaliasas.it



Scan the QR code to
discover more:
www.meditaliagroup.net/holdfinger/



If you need support about Holdfinger
contact us at **vendite.italia@meditaliasas.it**

Meditalia® | **endokey**

Meditalia Group S.r.l.
Via Saline, 75/A - 90149 Palermo (Italy)
C.F./P.I. 06932550822
Tel: +39 091 684 11 93
Email: info@meditaliagroup.net

www.meditaliagroup.net



Scopri di più su
www.meditaliagroup.net